

STOOL COLLECTION FOR OVA AND PARASITES (O&P)

COLLECT STOOL SPECIMEN IN A CLEAN, DRY CONTAINER. TRANSFER STOOL TO THE O&P KIT THAT WAS PROVIDED TO YOU BY YOUR PHYSICIAN OR LABORATORY STAFF USING THE FOLLOWING INSTRUCTIONS:

- 1. NAME AND BIRTHDATE SHOULD BE WRITTEN ON EACH VIAL.**
- 2. USING THE COLLECTION SPOON FROM THE KIT, PLACE A SMALL SCOOPFUL OF STOOL INTO EACH VIAL UNTIL THE FLUID LEVEL IN THE VIALS RISES TO THE RED ARROW INDICATED ON EACH VIAL.**
- 3. SECURELY REPLACE THE LID ON THE VIALS AND GENTLY SHAKE TO MIX THE CONTENTS.**
- 4. REFRIGERATE THE VIALS AND DELIVER TO THE LABORATORY AS SOON AS POSSIBLE.**

WARNING

THE LIQUID IN THE VIALS IS POISONOUS! PLEASE KEEP THEM OUT OF THE REACH OF CHILDREN. PLEASE CALL THE LABORATORY AT 717-273-6706 EXT 146 WITH ANY QUESTIONS.